



Curaplast® sensitive Wound Dressing Strips

Soft, supple nonwoven plaster for sensitive skin; first-aid dressing for hospital and doctor's surgery as well as for minor injuries at home and on the move.

The nonwoven material, which is slightly elastic widthwise, readily conforms to body contours and holds the dressing pad securely on the skin. A

micro-perforated separating film prevents adhesion in the wound area. The healing process is supported by dressing changes that are painless and gentle.



Curaplast® sensitive Wound Dressing Strips

In various sizes skin coloured nonwoven, 100% polyester, with continuous dressing pad, skin-friendly polyacrylate adhesive, free from colophonium and colophonium derivatives.

Curaplast® Sensitive Wound Dressing Strips			
for the practice and outpatient department	Size	REF	Shipping Units
5 m roll in cardboard dispenser	4 cm x 5 m	30 616	1 pc. / 36 pcs.
	6 cm x 5 m	30 617	1 pc. / 40 pcs.
	8 cm x 5 m	30 618	1 pc. / 24 pcs.
5 dressings, 10 cm each in cardboard box	6 cm x 0.5 m	17 076	1 pc. / 56 pcs.
10 dressings, 10 cm each	4 cm x 1 m	17 077	1 pc. / 47 pcs.
	6 cm x 1 m	17 078	1 pc. / 47 pcs.
	8 cm x 1 m	17 079	1 pc. / 36 pcs.
for household and travel			
Curaplast® Sensitive Finger Dressing			
for household and travel			
10 pcs. in cardboard box	2 x 18 cm	17 086	1 pack / 30 packs