

# Curaplast® **sensitive** Adhesive Wound Dressings



## Curaplast® **sensitive** Injection Plasters

Skin coloured precut adhesive plasters off the roll (cardboard dispenser); quick and easy to handle for blood tests, injections, etc.

## Curaplast® **sensitive** Precuts

Skin coloured, in different sizes, with central wound pad, each plaster sealed individually and hygienically.

<b>Curaplast® Sensitive Injection Plaster</b>			
in cardboard dispenser	Size	REF	Shipping Units
dispenser with 250 dressings	2 cm x 4 cm	30 625	1 pack / 16 packs
<b>Curaplast® Sensitive Precuts</b>			
<b>10 pcs., assorted sizes</b> in transparent pouch 6 dressings 4 dressings	1.9 cm x 7.2 cm 2.5 cm x 7.2 cm	17 082	1 pack / 340 packs
<b>20 pcs., assorted sizes</b> in cardboard box 10 dressings 4 dressings 4 dressings 2 dressings	1.9 cm x 7.2 cm 2.5 cm x 7.2 cm Ø 2.2 cm 3.8 cm x 3.8 cm	17 083	1 pack / 45 packs
<b>500 pcs.,</b> individually sealed in cardboard box	2.5 x 7.2 cm	30 631	1 pack / 38 packs
<b>Curaplast® Dressings, waterproof</b>		Skin coloured, precut adhesive plasters, 2-layer composite of polyurethane film and polyester nonwoven, skin-friendly polyacrylate adhesive.	
<b>10 pcs., assorted sizes</b> in transparent pouch 8 dressings 2 dressings	2.5 cm x 7.2 cm 1.9 cm x 7.2 cm	17 084	1 pack / 30 packs
<b>Curaplast® Kids</b>		Adhesive plaster with children's motifs; water and dirt repellent, sealed individually and hygienically.	
15 dressings in cardboard box	1.7 cm x 6 cm	30 630	1 pack / 48 packs
<b>Curaplast® Sensitive round</b>		Elastic, skin coloured nonwoven, with skin-friendly polyacrylate adhesive, free from colophonium and colophonium derivatives.	
in cardboard dispenser roll of 100 pcs.	Ø 2.3 cm	30 626	1 pack / 24 packs