



# Curapont® • Porofix® Wound Closure Strips

## Porofix® Wound Closure Strips

Windowed first-aid dressing, purpose-designed for adapting the wound edges in superficial stab wounds and incised wounds, or in lacerations; also suitable for relieving pressure on sutures.

In the windowed area the plaster has an aluminized film backing on the wound side, which prevents the plaster from sticking to the wound. Discharge is able to drain away, and the wound can be checked visually.

Unstretchable fabric with a high tensile strength, 100% viscose, skin coloured, coated with a skin-compatible synthetic rubber adhesive.



## Curapont® Wound Closure Strips

Plaster strips with high adhesive power, made from polyurethane reinforced with non-woven material, supple, tear-resistant and non-elastic; coated with hypoallergenic polyacrylate adhesive; for dependable, stitch-free wound closure; also for relieving pressure on sutures.

By adapting the wound edges cleanly and over a wide area, the resulting scar is perfect in cosmetic terms. The local anaesthetic

that is often required for a suture can be avoided; puncture channel infections, blood supply disorders in the suture area as well as the characteristic stitch marks are avoided.

Rapid and safe access can be gained at either end of the strip if necessary. This is particularly useful where half-length wound closure strips are sufficient; they are simply cut in half, complete with carrier paper, before being peeled off.



### Porofix® Wound Closure Strip

2 pcs. in a sealed bag

REF

Shipping Units

2.5 x 12.5 cm

30 100

10 pcs. / 1560 pcs.

### Curapont® Wound Closure Strips, sterile

sealed bags with

5 pcs.	3 x 75 mm	30 058	100 pcs. / 1000 pcs.
3 pcs.	6 x 75 mm	30 059	120 pcs. / 1200 pcs.
10 pcs.	6 x 100 mm	30 056	100 pcs. / 2000 pcs.
6 pcs.	12 x 100 mm	30 060	120 pcs. / 1200 pcs.
5 pcs.	12 x 140 mm	30 057	100 pcs. / 1000 pcs.