



Curapor® Island Dressing, sterile

Ready-for-use island dressing for initial sterile wound care and for changing dressings on wounds of all types – in the operating theatre, outpatient department and doctor's surgery. This dressing with its extremely high absorbing capacity has a polyester/polyethylene layer on the wound side; this allows discharge to diffuse into the absorptive part, accelerates wound

cleansing, yet does not stick to the wound. This has a positive influence on the healing process by ensuring that the wound is not disturbed when the dressing is changed. Reliable all-round closure provided by particularly adherent skin-friendly polyacrylate adhesive; the Curapor wound dressing can be removed from the skin painlessly and without leaving any residue.



Curapor® Island Dressing, sterile

Elastic nonwoven carrier material, white, 100% polyester, with skin-friendly polyacrylate adhesive (free from colophonium and colophonium derivatives); dressing pad with an absorbent layer made from viscose/polyester and a non-adhesive polyester/poly-ethylene wound contact layer, individually sealed.

Curapor® Island Dressing, sterile			
Size	Wound Pad	REF	Shipping Units
7 x 5 cm	4 x 2.5 cm	30 504	5 pcs. / 100 pcs.
10 x 8 cm	6 x 4 cm	30 505	5 pcs. / 100 pcs.
7 x 5 cm	4 x 2.5 cm	22 120	100 pcs. / 600 pcs.
10 x 8 cm	6 x 4 cm	22 121	50 pcs. / 300 pcs.
10 x 15 cm	5 x 11 cm	22 122	50 pcs. / 300 pcs.
10 x 20 cm	5 x 16 cm	22 123	50 pcs. / 300 pcs.
10 x 25 cm	5 x 20 cm	22 124	50 pcs. / 250 pcs.
10 x 30 cm	5 x 25 cm	22 125	50 pcs. / 250 pcs.
10 x 34 cm	5 x 29 cm	22 126	50 pcs. / 250 pcs.