




# Gazin® Gauze Balls

## **Gazin® Gauze Balls**

Made from DIN EN 14079 absorbent gauze 20 threads; optionally with or without X-ray contrast thread; for absorbing blood and discharge, for wound cleansing, also for disinfection of the skin.

## **Gazin® Dissecting Swabs**

Made from DIN EN 14079 absorbent gauze 24 threads; with X-ray contrast thread; firm construction for preparing layers of tissue and anatomical structures for surgical interventions.



<b>Gazin® Gauze Balls</b>		
	<b>non-sterile</b> without X-ray detectable thread in a pouch	
	REF	Shipping Units
walnut size	15 175	1000 pcs. / 5000 pcs.
plum size	15 176	1000 pcs. / 5000 pcs.
egg size	15 177	500 pcs. / 2500 pcs.
extra large	15 179	500 pcs. / 2000 pcs.
fist size	15 180	250 pcs. / 1000 pcs.
<b>non-sterile</b> with X-ray detectable thread in a pouch		
walnut size	15 183	1000 pcs. / 5000 pcs.
plum size	15 184	1000 pcs. / 5000 pcs.
egg size	15 185	500 pcs. / 2500 pcs.
extra large	15 186	500 pcs. / 2000 pcs.
fist size	15 223	250 pcs. / 1000 pcs.
<b>Gazin® Gauze Balls</b>		
	<b>sterile</b> in blister pack, packs with 2+3 or 2+2 pcs., without X-ray detectable thread	
plum size, 2+2 pcs.	14 961	100 pcs. (25x 4) / 1000 pcs. (250x 4)
plum size, 2+3 pcs.	14 962	125 pcs. (25x 5) / 1250 pcs. (250x 5)
egg size, 2+3 pcs.	14 968	100 pcs. (20x 5) / 1000 pcs. (200x 5)
extra large, 2+3 pcs.	14 969	100 pcs. (20x 5) / 1000 pcs. (200x 5)
dialysis, 2+2 pcs.	14 965	100 pcs. (25x 4) / 1000 pcs. (250x 4)
dialysis, 2+3 pcs.	14 966	125 pcs. (25x 5) / 1250 pcs. (250x 5)
with X-ray detectable thread		
plum size, 2+3 pcs.	14 967	125 pcs. (25x 5) / 1250 pcs. (250x 5)
<b>sterile</b> in blister pack, without X-ray detectable thread		
plum size, 3 pcs.	15 501	90 pcs. (30x 3) / 720 pcs. (240x 3)
plum size, 10 pcs.	15 500	300 pcs. (30x10) / 1500 pcs. (150x10)
plum size, 20 pcs.	15 499	400 pcs. (20x20) / 2000 pcs. (100x20)
egg size, 10 pcs.	15 498	300 pcs. (30x10) / 1500 pcs. (150x10)
<b>Gazin® Dissecting Swabs</b>		
	<b>non-sterile</b> with X-ray detectable thread in a pouch	
small	15 155	1000 pcs. / 3000 / 45000 pcs.
medium	15 156	1000 pcs. / 2000 / 30000 pcs.
large	15 157	500 pcs. / 1000 / 15000 pcs.