



Metalline® Non-adherent Dressing Kits and Dressing Cloths

Aluminized nonwoven material, sometimes combined with absorbent material, for the treatment of large area wounds and open wounds. The aluminized surface of Metalline dressing does not stick to the wound.

Metalline® Non-adherent Dressing Kit

For first aid and emergency care.

Combination of Metalline dressing 35 x 45 cm and gauze bandage 6 cm x 3 m.

Metalline® Non-adherent Bed Sheet

Aluminized nonwoven combined with absorbent material, for use as an underlay in the treatment of large-scale burns, chemical burns and abrasions.

Aluminized nonwoven material

combined with absorbent nonwoven material.

Metalline® Non-adherent Cloth

Aluminized nonwoven material off the roll, for use as an underlay or for covering large-scale wounds, for the cautious bedding of patients.

Aluminized nonwoven material off the roll. Metalline non-adherent cloth and any other absorbent layers must be sterilised before use.

Metalline® Non-adherent Dressing Cloth

For the care of large-scale injuries, burns, chemical burns and abrasions, in first aid and for emergency use.

Aluminized nonwoven combined with absorbent nonwoven material. The cloths are provided with ties to make them easier to handle during positioning and securing.



Metalline® Non-adherent Dressing Kit

| | REF | Shipping Units |
|--|--------|------------------|
| sterile , in foil pouch | | |
| consisting of: 1 dressing 35 cm x 45 cm 1 gauze bandage 6 cm x 3 m | 23 000 | 1 pc. / 150 pcs. |

Metalline® Non-adherent Bed Sheet

| | | |
|--------------------------------|--------|-----------------|
| sterile , in foil pouch | | |
| 2.5 m x 73 cm | 23 050 | 1 pc. / 20 pcs. |

Metalline® Non-adherent Cloth

| | | |
|---|--------|-----------------|
| non-sterile , rolled, in cardboard box | | |
| 10 m x 73 cm | 23 070 | 1 pc. / 10 pcs. |

Metalline® Non-adherent Dressing Cloth

| | | |
|--------------------------------|--------|------------------|
| sterile , in foil pouch | | |
| 40 x 60 cm | 23 012 | 1 pc. / 100 pcs. |
| 60 x 80 cm | 23 010 | 1 pc. / 100 pcs. |
| 80 x 120 cm | 23 011 | 1 pc. / 100 pcs. |