

Metalline® Wound Dressings

The aluminized wound pad of the Metalline dressings is a supple and soft, extremely fine non-woven material with good drainage. The smooth surface does not stick to the wound. This ensures that the wound is not disturbed, and the dressing can be changed without trauma, thus promoting the healing process. **Metalline Drainage Dressing:** the central slit and those radiating from the punch hole ensure a good fit around small calibre drains.



Metalline®

Multi-layered special dressing material.

Wound side: fine nonwoven veil, with vacuum-aluminized fibres.

Absorbent layer: high-bulk non-woven material.

Rear side: needle-punched non-woven material. Layers joined through needle punching and impervious dispersion adhesive, resulting in good absorptive characteristics.



Metalline® Dressings

For all types of wounds – in traumatology, surgery, dermatology and phlebology, especially for extensive wounds, abrasions, burns; also for first aid use.

Metalline® Drainage Dressings

For care of wounds when small calibre drains or infusion catheters are involved, as well as for wire or nail extensions.

Hole diameter:
from 4.5 mm - 11 mm
for Charrière 9 - 30

Metalline® Tracheo Dressings

For care with tracheostomies or large-diameter drains.

Hole diameter:
from 12 mm - 19 mm
for Charrière 31 - 50

Metalline® Wound Dressings			
	Size	REF	Shipping Units
sterile, individually sealed	8 x 10 cm	23 083	10 pcs. / 560 pcs.
	10 x 12 cm	23 084	10 pcs. / 300 pcs.
sterile, individually sealed; for OR, outpatient department and practice, in dispenser box	8 x 10 cm	23 081	50 pcs. / 600 pcs.
	10 x 12 cm	23 082	50 pcs. / 300 pcs.
	10 x 20 cm	23 087	50 pcs. / 300 pcs.
non-sterile, rolled, in foil pouch	5 m x 10 cm	23 080	1 pc. / 20 pcs.
Metalline® Drainage Dressing			
sterile, individually sealed, in dispenser box	6 x 7 cm	23 093	50 pcs. / 600 pcs.
Metalline® Tracheo Dressing			
sterile, individually sealed, in dispenser box	8 x 9 cm	23 094	50 pcs. / 600 pcs.