

Suprasorb® F

Film Wound Dressing

Suprasorb® F Protects and holds (in)visibly!

Suprasorb F offers transparent help for protection and dressing retention. With the breathability of a semi-permeable polyurethane film.

Suprasorb F enables a moist wound climate conducive to wound healing by means of gas

and heat exchange, repels fluids and protects the wound against colonisation by germs from outside.


Indications: lightly exuding, superficial, non-infected wounds in the epithelialisation phase, e.g.: abrasions, dermabrasion, ulcers, decubitus, diabetic foot syndrome (DFS), second degree burns.



Similar to skin and extremely fine.

The Suprasorb F microporous film wound dressing offers a variety of "clear" benefits: visual wound inspection, very easy application and use, effective protection from external influences, allows showering and washing, flexible, elastic and tear-resistant.

Suprasorb F is available as a sterile film wound dressing and also as a non-sterile film dressing off the roll.



- 1 Transparent, semi-permeable polyurethane film.
- 2 Repels bacteria and fluids.
- 3 Surplus water vapour evaporates.

Interval between dressing changes

Depending on the degree of exudation the dressing should be changed when there is an accumulation of exudate over the wound size.

Suprasorb® F Film Wound Dressing, sterile, individually sealed

Adhesive Surface	REF	Shipping Units
5 x 7 cm	20 460	10 pcs. / 380 pcs.
5 x 7 cm	20 461	100 pcs. / 1000 pcs.
10 x 12 cm	20 462	10 pcs. / 100 pcs.
10 x 12 cm	20 463	50 pcs. / 200 pcs.
10 x 25 cm	20 464	10 pcs. / 140 pcs.
15 x 20 cm	20 465	10 pcs. / 400 pcs.
20 x 30 cm	20 466	10 pcs. / 180 pcs.

Suprasorb® F Film Wound Dressing, off the roll, non-sterile

10 cm x 1 m	20 467	1 pc. / 16 pcs.
10 cm x 10 m	20 468	1 pc. / 16 pcs.
15 cm x 10 m	20 469	1 pc. / 10 pcs.