



# Vliwazell®

## Absorbent Wound Dressing

Fully-enclosed, supple dressing with very high absorption capacity; for the care of wounds with heavy discharge. Vliwazell dressings are wound-friendly

through their particularly soft covering of nonwoven material. Wound exudate is quickly channelled into the absorptive cellulose layer, whose capacity is

many times greater than that of a gauze dressing. A moisture inhibiting layer of nonwoven material on the back protects the clothing.



Vliwazell®	Size	REF	Shipping Units
<b>non-sterile,</b> 25 pcs. in plastic bag	10 x 10 cm	30 430	25 pcs. / 900 pcs.
	10 x 20 cm	30 431	25 pcs. / 450 pcs.
	15 x 20 cm	30 432	25 pcs. / 200 pcs.
	15 x 25 cm	30 433	25 pcs. / 200 pcs.
	20 x 20 cm	30 434	25 pcs. / 200 pcs.
	20 x 25 cm	30 435	25 pcs. / 150 pcs.
	20 x 40 cm	30 436	25 pcs. / 200 pcs.
<b>sterile,</b> individually sealed	10 x 10 cm	30 450	60 pcs. / 360 pcs.
	10 x 20 cm	30 451	30 pcs. / 300 pcs.
	15 x 20 cm	30 452	30 pcs. / 180 pcs.
	15 x 25 cm	30 453	30 pcs. / 240 pcs.
	20 x 20 cm	30 454	30 pcs. / 120 pcs.
	20 x 25 cm	30 455	30 pcs. / 150 pcs.
	20 x 40 cm	30 456	30 pcs. / 180 pcs.