

Soft, economical universal dressing, made from a combination of a fine nonwoven layer and bleached cellulose; for wounds with heavy discharge, as an additional ab-

sorbent layer, for ointment dressings and moist dressings; also for padding.

The wound-friendly nonwoven material and the absorbent lay-

er are joined together without any creases. Zemuko dressings can thus be cut off the roll with well-defined edges and no fluff.

Zemuko®

Cellulose/nonwoven combination. Wound side: nonwoven, 67% viscose, 33% polyester; absorbent layer: bleached cellulose (in accordance with German Pharmacopoeia DAB); the two layers are joined together by embossing; sterilised, or sterilisable (A 121°C).

Zemuko® seamless

Alternative universal dressing; made from bleached cellulose DAB with gauze **covering on all sides**.



Zemuko®			
	Size	REF	Shipping Units
rolled, in plastic dispenser	10 cm x 1 m	23 512	1 pc. / 60 pcs.
	10 cm x 2 m	23 515	1 pc. / 30 pcs.
rolled, in a pouch for the practice and outpatient department	10 cm x 10 m	23 502	1 pc. / 20 pcs.
	15 cm x 10 m	23 503	1 pc. / 16 pcs.
	28 cm x 10 m	23 602	1 pc. / 10 pcs.
Zemuko® seamless			
rolled, in a pouch	10 cm x 10 m	23 702	1 pc. / 20 pcs.
	15 cm x 10 m	23 703	1 pc. / 16 pcs.
Zemuko® Wound Dressings sterile, individually sealed			
	10 x 10 cm	23 802	10 pcs. / 200 pcs.
Dispenser		Polystyrol, transparent, suitable for rolled dressing materials, in surgeries and outpatient departments, e.g. Gazin, Zemuko, Zelletten, etc. max. width of contents 10 cm	
empty, max. width of contents: 10 cm		23 900	1 pc.