

# Lomatuell® H Paraffin Gauze Dressings • Paraffin Ribbon Gauze

## Lomatuell® H

Ready-for-use tulle, impregnated with a hydrophobic neutral ointment base; does not stick to the wound; Lomatuell is a soothing dressing for surface wounds and burns; for radiation injuries and leg ulcers; for

dressings following skin grafts and dermatological indications; also postoperatively in plastic surgery.

The wide-meshed tulle promotes the drainage of discharge and ensures that the wound is well ventilated. The neutral oint-


ment base provides the opportunity to apply additional, wound-specific medication.

### Carrier material:

Wide-meshed tulle, 100% cotton;

### Ointment base:

White paraffin



## Paraffin Ribbon Gauze

Paraffin-impregnated ribbon gauze for packing cavity wounds and fistulae, also in the fields of gynaecology and proctology; following surgery and for application in the fields of ear, nose and throat therapy and orthodontic surgery.

### Carrier material:

24-thread ribbon gauze with selvedges, 100% cotton;

### Ointment base:

White paraffin

<b>Lomatuell® H Paraffin Gauze Dressing, sterile</b>			
	Size	REF	Shipping Units
individually sealed	5 x 5 cm	23 314	10 pcs. / 360 pcs.
	10 x 10 cm	23 315	10 pcs. / 120 pcs.
	10 x 20 cm	23 316	10 pcs. / 120 pcs.
	10 x 30 cm	23 317	10 pcs. / 80 pcs.
	10 x 10 cm	23 318	50 pcs. / 450 pcs.
individually in plastic box	10 cm x 7 m	93 903	1 pc. / 72 pcs.
<b>Paraffin Ribbon Gauze</b>			
individually in foil pouch	2 cm x 1 m	13 514	10 pcs. / 300 pcs.
individually in plastic jar	2 cm x 10 m	13 515	1 pc. / 96 pcs.