

# Opraclean® Dressing Materials

## Gazin® Ribbon Gauze

### Opraclean® Dressing Materials

Opraclean supports wound cleansing by absorbing wound exudate, cellular debris and bacteria.

Opraclean is well suited for packing even deep wound pockets and cavities. At the same time, the dressings which are impregnated with iodoform have an antimicrobial effect. When it comes into contact with wound exudates, iodoform dissolves and continuously releases small quantities of iodine from the dressing. Iodine is known for its broad effective spectrum against a multitude of gram-negative, gram-positive, aerobic and anaerobic bacteria and is also effective against many fungi and viruses.

For antiseptic packing and drainage after surgical debridement in infected wounds and wounds with a risk of infection, for example abscesses, fistulae, ulcers, and in body cavities. For acute and chronic wounds.

Handy dispensers enable easy handling in everyday clinical and

general practice and satisfy hygienic requirements at the same time.

Carrier material 100% gauze, impregnated with Iodoform.

### Gazin® Ribbon Gauze

Ribbon gauze  
DIN 61 635-TB (carrier material)



<b>Opraclean® Ribbon Gauze with Iodoform, off the roll</b>			
	Size	REF	Shipping Units
in a dispenser, in a special pouch	1 cm x 5 m	13 580	1 pc. / 216 pcs.
	2 cm x 5 m	13 581	1 pc. / 120 pcs.
	5 cm x 5 m	13 582	1 pc. / 60 pcs.
<b>Opraclean® Tubular Gauze with Iodoform</b>			
in a dispenser, in a special pouch	8 cm x 5 m	13 583	1 pc. / 30 pcs.
<b>Opraclean® Absorbent Gauze with Iodoform, 4-ply</b>			
in a dispenser, in a special pouch	6 cm x 5 m	13 584	1 pc. / 30 pcs.
	8 cm x 5 m	13 585	1 pc. / 30 pcs.
<b>Gazin® Ribbon Gauze</b>			
sterile, rolled, in transparent dispenser	1 cm x 5 m	13 470	1 pc. / 216 pcs.
	2 cm x 5 m	13 471	1 pc. / 120 pcs.
	5 cm x 5 m	13 473	1 pc. / 60 pcs.