

Vliwaktiv®

Vliwaktiv® Ag

Vliwaktiv® Activated Charcoal Dressing

Ideal combination of a highly absorbent pad with a wound dressing that absorbs both odour and bacteria, for the treatment of malodorous wounds such as leg ulcers, decubital ulcers, anal fistulae, etc.

A skin-friendly, breathable non-woven cover speeds up the passage of exudate into the absorbent core. Bacteria and odours are effectively bound in Vliwaktiv. The clothing protection prevents the seeping of wound exudate.

Vliwaktiv® Ag Activated Charcoal Dressing with Silver

Activated charcoal dressing impregnated with silver, in a high-quality nonwoven cover that does not stick to the wound; for the treatment of malodorous infected and critically colonised wounds, such as leg ulcers, decubital ulcers, exulcerative tumors, etc. Presentations: rope and absorbent dressings.

Vliwaktiv Ag combines the antimicrobial action of silver with the benefits of activated charcoal which adsorbs odours and binds bacteria.

When Vliwaktiv Ag makes contact with exudate, it releases silver ions and destroys the retained microorganisms. Thus, microbial counts are effectively reduced without undue stress to the wound.

Unpleasant wound odours which are often associated with infected wounds are effectively adsorbed by the silver/activated

charcoal fabric, which enhances the well-being of patients and carers alike. Vliwaktiv Ag also effectively binds endotoxins which may impair the wound healing process.

Deep wounds and wound cavities can be treated easily with the flexible, soft Vliwaktiv activated charcoal rope. For heavily exuding superficial wounds, the Vliwaktiv Ag activated charcoal dressing is made with a highly absorbent core and integrated clothing protection which prevents the seeping of exudate.



Vliwaktiv® sterile			
	Size	REF	Shipping Units
individually sealed	10 x 10 cm	20 254	20 pcs. / 400 pcs.
	10 x 20 cm	20 255	20 pcs. / 200 pcs.
Vliwaktiv® Ag Activated Charcoal Rope with Silver, sterile			
individually sealed	6.5 x 10 cm	20 700	10 pcs. / 300 pcs.
	10 x 10 cm	20 701	10 pcs. / 300 pcs.
	10 x 15 cm	20 702	10 pcs. / 340 pcs.
Vliwaktiv® Ag Absorbent Activated Charcoal Dressing with Silver, sterile			
individually sealed	10 x 10 cm	20 705	10 pcs. / 150 pcs.
	10 x 20 cm	20 706	10 pcs. / 110 pcs.
	20 x 20 cm	20 707	10 pcs. / 130 pcs.