

# Cannula Retention

Special dressings available in various materials and forms for reliable and fast retention of indwelling cannulas, butterfly cannulas and drains.

## Curapor® i.v.

Nonwoven retention dressing with punched slit for smooth and secure positioning of indwelling cannulas; with a separate cushion for swabbing/padding.

Fine, elastic nonwoven, white, 100% polyester, skin-friendly polyacrylate adhesive, indifferent to radiation.



## Porofix® Cannula Retention

Skin coloured adhesive tape for retaining various types of infusion catheters and their tubes. Non-elastic fabric (100% viscose) with punched central hole/slit; Skin-friendly synthetic rubber adhesive, highly air-permeable.



## Suprasorb® F Film Wound Dressing

This transparent, extremely fine polyurethane film is elastic and supple like the skin itself; coated with a skin-friendly polyacrylate adhesive.



## Curagard®

Transparent, breathable membrane with foam edging and an elastic strip for the retention of peripheral and central venous catheters. The membrane forms a barrier against microorganisms and fluids and allows visual inspection of the puncture point at any time.

The dressing conforms well to body contours, adheres reliably to the skin and protects the area around the puncture point against germs and bacteria. Suprasorb F is waterproof, but – thanks to its permeability to water vapour – does not impede skin functions.

## Curafix® i.v. control

Nonwoven/film laminate with a transparent film over the puncture site for monitoring. For the retention of peripheral venous and arterial indwelling cannulas. Easy application, good fit providing reliable protection against infection, can be wiped, resistant to disinfectants, latex-free. High water vapour permeability, with a highly skin-friendly polyacrylate adhesive.

With an additional sterile pad.



<b>Curapor® i.v. sterile</b>			
	Size	REF	Shipping Units
individually sealed	6.5 x 8.8 cm	30 520	50 pcs. / 400 pcs.
<b>Curagard® sterile, individually sealed</b>			
<b>Curagard® SP</b> omega shaped, for adults	7.5 x 6.5 cm	30 117	100 pcs. / 500 pcs.
<b>Curagard® JR</b> omega shaped, for children	4 x 6.5 cm	30 116	100 pcs. / 500 pcs.
<b>Curagard® AP</b> rectangular	7 x 6 cm	30 115	100 pcs. / 500 pcs.
<b>Curagard® TPN</b> rectangular	13 x 10 cm	30 118	25 pcs. / 150 pcs.
<b>Curafix® i.v. control sterile</b>			
individually sealed	9 x 6 cm	30 200	50 pcs. / 400 pcs.
<b>Curafix® i.v. classic Cannula Retention</b>			
	2.5 x 12.5 cm	30 102	20 pcs. / 2000 pcs.
<b>Suprasorb® F Film Wound Dressing sterile</b>			
individually sealed	5 x 7 cm 5 x 7 cm	20 460 20 461	10 pcs. / 380 pcs. 100 pcs. / 1000 pcs.

