



Crimpelast

Raucolast Conforming Bandages

Crimpelast

Skin-friendly, highly absorbent elastic gauze bandage with good extensibility; allows easy and time saving application.

Even when joints are moved a lot, the bandage remains in place without slipping. Crimpelast is suitable for all types of dressing retention.



Crimpelast		White, 100% cotton, lengthwise elastic, extensibility approx. 40%		
	Width	REF	Shipping Units	
4 m (stretched length) 20 per box	6 cm	10 101	1 box / 10 boxes (200 pcs.)	
	8 cm	10 102	1 box / 10 boxes (200 pcs.)	
	10 cm	10 099	1 box / 8 boxes (160 pcs.)	
	12 cm	10 100	1 box / 6 boxes (120 pcs.)	
	15 cm	13 997	1 box / 5 boxes (100 pcs.)	
Raucolast		With selvages, white, 52% polyamide, 48% viscose, elastic, extensibility approx. 110%, washable, sterilisable (A 121°C)		
4 m individually sealed	6 cm	11 495	20 pcs. / 200 pcs.	
	8 cm	11 498	20 pcs. / 200 pcs.	
	10 cm	11 484	20 pcs. / 160 pcs.	
	12 cm	11 487	20 pcs. / 120 pcs.	
20 per box	4 cm	11 493	1 box / 16 boxes (320 pcs.)	
	6 cm	11 496	1 box / 10 boxes (200 pcs.)	
	8 cm	11 499	1 box / 10 boxes (200 pcs.)	
	10 cm	11 485	1 box / 8 boxes (160 pcs.)	
	12 cm	11 488	1 box / 6 boxes (120 pcs.)	