



Geka® Gauze Bandage

Geka gauze bandages are used mainly for dressing retention purposes and, in connection with splints, for the immobilisation of injured extremities; to close split casts or as a final layer over zinc paste bandages (Unna's Boot).



White, with selvages, 20 threads, 100% viscose, sterilisable (A 121 °C).

Geka® Gauze Bandage			
	Size	REF	Shipping Units
4 m unwrapped 20 per box	4 cm	19 160	20 pcs. / 880 pcs.
	6 cm	19 161	20 pcs. / 600 pcs.
	8 cm	19 162	20 pcs. / 480 pcs.
	10 cm	19 163	20 pcs. / 400 pcs.
	12 cm	19 164	20 pcs. / 320 pcs.