



Lenkelast®

Elastra Ideal medium stretch

Firm-grip, all purpose bandages with durable elasticity and with a wide range of applications: everyday retention bandages with compressive wound care, or splint dressings, for support

and pressure relief; also for light compression.

The extensibility of these bandages makes them conform well to all body contours, and the dressing remains in posi-

tion. The fine fabric is air-permeable and highly skin-friendly.

Lenkelast®

White, 64% cotton, 35% polyamide, 1% elastane, medium stretch, extensibility approx. 130%.

Elastra Ideal

White, 58% cotton, 41% polyamide, 1% elastane, medium stretch, extensibility approx. 135%, washable, sterilisable (A 121°C).



Lenkelast®			
	Width	REF	Shipping Units
5 m (stretched length) individually sealed	6 cm	19 460	10 pcs. / 250 pcs.
	8 cm	19 461	10 pcs. / 200 pcs.
	10 cm	19 462	10 pcs. / 240 pcs.
	12 cm	19 463	10 pcs. / 200 pcs.
	15 cm	19 464	10 pcs. / 160 pcs.
	20 cm	19 465	10 pcs. / 100 pcs.
5 m (stretched length) unwrapped 10 per box	6 cm	19 490	10 pcs. / 250 pcs.
	8 cm	19 491	10 pcs. / 200 pcs.
	10 cm	19 492	10 pcs. / 240 pcs.
	12 cm	19 493	10 pcs. / 200 pcs.
	15 cm	19 494	10 pcs. / 160 pcs.
	20 cm	19 495	10 pcs. / 100 pcs.
Elastra Ideal			
5 m (stretched length) unwrapped 10 per box	6 cm	10 122	10 pcs. / 160 pcs.
	8 cm	10 124	10 pcs. / 120 pcs.
	10 cm	10 108	10 pcs. / 120 pcs.
	12 cm	10 110	10 pcs. / 100 pcs.
	15 cm	10 116	10 pcs. / 80 pcs.
	20 cm	10 118	10 pcs. / 60 pcs.