

Mollelast®

Conforming Bandage

Firm-grip, highly elastic bandage for all types of dressing retention; the special knitting technique makes Mollelast widthwise stable and air-permeable. Mollelast is easy to handle and comfortable on the

skin; it conforms well to body contours without creasing and secures dressings reliably without cutting in, especially on cone-shaped body parts and in areas close to joints.



Knitted bandage, white, 56% viscose, 44% polyamide, elastic, extensibility approx. 100%, sterilisable (A 121°C).



Mollelast®			
	Size	REF	Shipping Units
4 m (stretched length) individually sealed	4 cm	14 410	20 pcs. / 520 pcs.
	6 cm	14 411	20 pcs. / 360 pcs.
	8 cm	14 412	20 pcs. / 500 pcs.
	10 cm	14 413	20 pcs. / 360 pcs.
	12 cm	14 414	20 pcs. / 320 pcs.
20 per box	4 cm	19 410	20 pcs. / 520 pcs.
	6 cm	19 411	20 pcs. / 360 pcs.
	8 cm	19 412	20 pcs. / 500 pcs.
	10 cm	19 413	20 pcs. / 360 pcs.
	12 cm	19 414	20 pcs. / 320 pcs.
100 per box	6 cm	19 427	100 pcs. / 400 pcs.
	8 cm	19 428	100 pcs. / 400 pcs.
	10 cm	19 429	100 pcs. / 300 pcs.