

tg® Ready-to-wear Dressings

These ready-to-wear dressings are ideal for covering body parts affected by dermatological disorders, e. g. contact dermatitis, photodermatitis and medically treated skin.

For such indications, ready-to-wear dressings are quick and easy to apply. The material used was selected for its excellent hypoallergenic characteristics.

tg® Glove

Soft, supple hand/finger dressing shaped for a comfortable fit.



Also for preventing superinfections.

100% cotton, washable, sterilisable (A 121°C).

tg® Mitten

Thumbless hand dressing with integrated tie strings for securing at the wrist. As protective covering for dermatological disorders, and to prevent children with itchy skin conditions from scratching themselves.

100% cotton, washable, sterilisable (A 121°C).



tg® Fingerstall

Ready-to-wear finger dressing (size 2) for securing dressings and for protective covering in the case of discharging, acutely inflammatory or chronic skin disorders.

67% cotton, 33% viscose, sterilisable (A 121°C).

tg® Head Dressing

Ready-to-wear head dressing (size 9)

67% cotton, 33% viscose, sterilisable (A 121°C).

tg® Hand-/Foot Dressing

Ready-to-wear hand-/foot dressing (size 7)

67% cotton, 33% viscose, sterilisable (A 121°C).

tg® Glove		REF	Shipping Units
each pair individually in plastic bag			
for children		24 749	1 pair / 100 pairs
for adults	Gr. 6-7 (small)	24 750	1 pair / 50 pairs
	Gr. 7½-8½ (medium)	24 751	1 pair / 50 pairs
	Gr. 9-10 (large)	24 752	1 pair / 50 pairs
tg® Mitten			
each pair individually small (infants)		24 746	1 pair / 100 pairs
in plastic bag large (small children)		24 747	1 pair / 100 pairs
tg® Fingerstall			
individually in pouch		24 080	1 pc. / 250 pcs.
50 pcs. in cardboard box		24 081	50 pcs. / 750 pcs.
tg® Head Dressing			
individually in pouch		24 082	1 pc. / 160 pcs.
10 pcs. in cardboard box		24 083	10 pcs. / 200 pcs.
tg® Hand-/Foot Dressing			
individually in pouch		24 084	1 pc. / 200 pcs.
10 pcs. in cardboard box		24 085	10 pcs. / 200 pcs.