

tg[®] Tubular Bandage

Practical and time-saving (tubular) bandaging technique – derived from the characteristics of seamlessly knitted material – for securing dressings and for light compression; even allows suspension and extension (with traction ring and traction ring belts).

The close-meshed, particularly skin-friendly material can be stretched to about four times its original width, and by pulling lengthwise can be narrowed down again; it is ideal for use next to the skin.

tg tubular bandage does not slip, but remains firmly in position without cutting in. The dressing remains comfortably smooth even on body parts that are difficult to dress. tg tubular bandage is a match for any dressing in any location from head to toe.

tg[®] Tubular Bandage

Seamlessly knitted tube made from fine yarns, white, 67% cotton, bleached, 33% viscose, particularly skin-friendly, economical in use, can be sterilised if necessary (A 121°C).

tg tubular bandage is also used for dermatological indications, over ointment dressings or zinc paste bandages or as undercast protection for P.o.P. or synthetic casts.



tg [®] Tubular Bandage			REF	Shipping Units
5 m roll white, individually in cardboard dispenser	size 1 (1.4 cm)	for finger and toe	24 020	1 roll / 80 rolls
	size 2 (2.3 cm)	for fingers and toes	24 021	1 roll / 50 rolls
	size 3 (3 cm)	for several fingers, child's arm	24 022	1 roll / 40 rolls
	size 5 (5.5 cm)	for arm, lower leg, child's leg	24 023	1 roll / 40 rolls
	size 6 (6.5 cm)	for arm, leg (medium)	24 024	1 roll / 40 rolls
	size 7 (7 cm)	for arm (large), leg, child's head	24 025	1 roll / 40 rolls
	20 m roll	size 1 (1.4 cm)	for finger and toe	24 000
size 2 (2.3 cm)		for fingers and toes	24 001	1 roll / 20 rolls
size 3 (3 cm)		for several fingers, child's arm	24 002	1 roll / 10 rolls
size 5 (5.5 cm)		for arm, lower leg, child's leg	24 003	1 roll / 10 rolls
size 6 (6.5 cm)		for arm, leg (medium)	24 004	1 roll / 10 rolls
size 7 (7 cm)		for arm (large), leg, child's head	24 005	1 roll / 18 rolls
size 9 (8.5 cm)		for thigh (large), head, armpit	24 006	1 roll / 12 rolls
size 12 (12 cm)	for thigh (extra large), head (large), child's trunk	24 007	1 roll / 12 rolls	
10 m roll	size K1 (16 cm)	for trunk (small) < M	24 030	1 roll / 10 rolls
	size K2 (21 cm)	for trunk (large) > L	24 031	1 roll / 10 rolls