



## tg® Vest

Short-sleeved vest, soft, supple and hypoallergenic, as a practical ready-to-wear dressing for covering extensive skin disorders; also as a time-saving under-cast protection for P.o.P. and synthetic casts.

The elastic knitted material has widthwise extensibility allowing for intermediate sizes.

Ribbed material with widthwise extensibility, raw white, 100% cotton, not degreased or optically brightened, particularly hypoallergenic, washable, sterilisable (A 121°C).



### tg® Vest

	REF	Shipping Units
individually in cardboard box		
<b>small</b> (ladies' size 36 - 44, men 42 - 50)	24 044	1 pc. / 14 pcs.
<b>large</b> (ladies' up to size 46, men up to size 56)	24 045	1 pc. / 9 pcs.