

tg[®] fix Tubular Net Bandage

Wide-mesh, highly elastic tubular net bandage for fast and easy retention of all types of wound dressings, also for the application of occlusion films. tg fix conforms well to all body contours without creasing, it does not tighten and is hardly felt. The material is ladder-

proof; the dressing can be cut open at any point without the cuts tearing further. For wound inspection and dressing change it is usually sufficient to partially lift the tubular net or to pull it back. For this reason it is also well suited for aftercare at home.

The net-type, knitted tube is made from fine polyamide-covered rubber threads and crimped polyamide yarn; white, free from optical brighteners; 71% polyamide, 29% elastodiene, bi-elastic, can be washed and sterilised (A 121°C).



tg[®] fix Tubular Net Bandage

white

4 m (stretched length)

folded, in a cardboard box

Size

REF

Shipping Units

for finger

for small extremities (several fingers together/hand/foot)

for large extremities (small head/arm/leg)

for large head, small trunk

for large trunk, hip, armpit

A

24 240

1 pc. / 100 pcs.

B

24 241

1 pc. / 100 pcs.

C

24 242

1 pc. / 98 pcs.

D

24 243

1 pc. / 58 pcs.

E

24 244

1 pc. / 28 pcs.

25 m (stretched length)

folded, in a cardboard dispenser

for finger

for small extremities (several fingers together/hand/foot)

for large extremities (small head/arm/leg)

for large head, small trunk

for large trunk, hip, armpit

A

24 250

1 pc. / 66 pcs.

B

24 251

1 pc. / 20 pcs.

C

24 252

1 pc. / 16 pcs.

D

24 253

1 pc. / 16 pcs.

E

24 254

1 pc. / 12 pcs.