

Dauerbinde® • Dauerbinde® biforte Perfekta® long stretch

Dauerbinde® K and F

Long stretch bandages with strong compression, with a characteristic firm-grip fabric texture in Dauerbinde K and a finer weave in Dauerbinde F.

Due to their long stretch properties, the compressive power of these bandages can be easily adjusted. They are used as support and relief bandages on ligaments and joints.

Because of the high resting pressure the bandage must be removed at night or when the patient is in a resting position.

The Dauerbinde can be fitted easily even without much experience in bandaging techniques; the bandages conform well to body contours, have little tendency to bulk up, are durable and easy-care.

Skin coloured, durable elasticity, long stretch, extensibility approx. 180%, washable at 95°C.



Dauerbinde® biforte

Two-way stretch elastic bandage. The bandage is indicated in sports injuries, in supporting and relieving sprains and contusions, for after-care in traumatology, and for strong and permanent compression in phlebology.

Skin coloured, bi-elastic long stretch, 144% in length, 84% in width, washable at 95°C.



Perfekta® Super (extra fine)

Long stretch bandage with durable elasticity, with low power of return; particularly well tolerated by the skin due to its high cotton content.

Support and relief bandage for prophylactic relief of capsules and ligaments, prevention of haematomas and for contusions, distortions, dislocations and tendon injuries.

Perfekta Super can also be used as a compressive bandage.

Skin coloured, resistant to ointments, long stretch, extensibility approx. 160%, washable at 95°C.

Perfekta® Kräftig (strong)

Long stretch bandage with durable elasticity and strong power of return; particularly well tolerated by the skin due to its high cotton content.

Compression bandages for varicosis, oedemas, angiodysplasia, sclerotherapy, following surgical interventions on veins and for prevention of thrombosis. Support and relief bandages for phlebitis or existing injuries.

Perfekta® Fein (fine)

Same as Perfekta Kräftig, but with a ribbed surface and medium power of return.

Skin coloured, resistant to ointments, long stretch, extensibility approx. 170%, washable at 95°C.



Dauerbinde® K		86% cotton, 7% polyamide, 7% elastane		
	Width	REF	Shipping Units	
7 m (stretched length) with bandage clips individually in cardboard box	6 cm	22 000	1 pc. / 100 pcs.	
	8 cm	22 001	1 pc. / 100 pcs.	
	10 cm	22 002	1 pc. / 100 pcs.	
	12 cm	22 003	1 pc. / 50 pcs.	
	20 cm	22 004	1 pc. / 20 pcs.	
14 m (stretched length) with bandage clips individually sealed	10 cm	22 007	1 pc. / 25 pcs.	
	12 cm	22 008	1 pc. / 25 pcs.	
7 m (stretched length) with bandage clips unwrapped 10 per box	6 cm	22 020	10 pcs. / 200 pcs.	
	8 cm	22 021	10 pcs. / 180 pcs.	
	10 cm	22 022	10 pcs. / 120 pcs.	
	12 cm	22 023	10 pcs. / 100 pcs.	
	20 cm	22 024	10 pcs. / 60 pcs.	
Dauerbinde® F		82% cotton, 13% polyamide, 5% elastane		
7 m (stretched length) with bandage clips individually in cardboard box	6 cm	22 010	1 pc. / 100 pcs.	
	8 cm	22 011	1 pc. / 100 pcs.	
	10 cm	22 012	1 pc. / 100 pcs.	
	12 cm	22 013	1 pc. / 50 pcs.	
7 m (stretched length) with bandage clips unwrapped 10 per box	6 cm	22 025	10 pcs. / 160 pcs.	
	8 cm	22 026	10 pcs. / 80 pcs.	
	10 cm	22 027	10 pcs. / 120 pcs.	
	12 cm	22 028	10 pcs. / 80 pcs.	
Dauerbinde® biforte		69% viscose, 24% cotton, 7% elastodiene		
7 m (stretched length) individually sealed	8 cm	16 980	1 pc. / 40 pcs.	
	10 cm	16 981	1 pc. / 30 pcs.	
	12 cm	16 982	1 pc. / 28 pcs.	
Perfekta® Super (extra fine)		96% cotton, 4% elastane		
7 m (stretched length) unwrapped 8 per box	8 cm	20 136	8 pcs. / 80 pcs.	
	10 cm	20 131	8 pcs. / 80 pcs.	
	12 cm	20 133	8 pcs. / 80 pcs.	
Perfekta® Kräftig (strong)		94% cotton, 6% elastane		
7 m (stretched length) unwrapped 8 per box	8 cm	20 129	8 pcs. / 80 pcs.	
	10 cm	20 109	8 pcs. / 80 pcs.	
	12 cm	20 117	8 pcs. / 80 pcs.	
Perfekta® Fein (fine)		94% cotton, 6% elastane		
7 m (stretched length) individually in cardboard box	8 cm	20 103	1 pc. / 50 pcs.	
	10 cm	20 099	1 pc. / 50 pcs.	
	12 cm	20 100	1 pc. / 50 pcs.	
10 m (stretched length)	10 cm	20 098	1 pc. / 50 pcs.	
	12 cm	93 185	1 pc. / 50 pcs.	