

Dauerlastic® TPS

Thrombo-prophylactic Stockings

Seamless, circular knitted stocking for thrombo-prophylaxis.

The stocking is widthwise elastic and shaped anatomically in such a way that the compression applied to the leg and the veins in the leg falls gradually from 18 mmHg at the ankle to a

maximum of 8 mmHg at the thigh. The foot and ankle are firmly enclosed.

Dauerlastic TPS are provided with an inspection opening below the toes. Dauerlastic TPS stockings assist the blood to flow back through the veins in

immobilised patients before, during and after operations, and in obstetrics.

Seamless, circular knitted, white, 80% polyamide, 20% elastane; flat garter strip with silicone rubber coating, can be boiled and sterilised (A 121°C).



The compression values are unaffected even by repeated washing.

The stocking sizes fit approx. 90% of patients; if the given values fall below or above the specified dimensions, compression bandages (e. g. Dauerbinde K or F) should be used instead.

Contraindications: arterial disorders which act haemodynamically; severe cardiac insufficiency; serious oedemas in the leg; dimensions of the thigh which exceed the maximum dimensions of the stockings; care should be exercised with venous ligatures.

Dauerlastic® TPS

Thigh length, one pair per pack

Size	Colour code foot	Thigh length to gluteal fold	Calf circumference	Thigh circumference (max.)	REF	Shipping Units
A1	1 blue stripe	up to 75 cm	up to 34.5 cm	up to 65 cm	23 570	1 pack / 26 packs
A2	2 blue stripes	75-85 cm	up to 34.5 cm	up to 65 cm	23 571	1 pack / 26 packs
A3	3 blue stripes	over 85 cm	up to 34.5 cm	up to 65 cm	23 572	1 pack / 26 packs
B1	1 green stripe	up to 75 cm	34.5-40.5 cm	up to 70 cm	23 573	1 pack / 26 packs
B2	2 green stripes	75-85 cm	34.5-40.5 cm	up to 70 cm	23 574	1 pack / 26 packs
B3	3 green stripes	over 85 cm	34.5-40.5 cm	up to 70 cm	23 575	1 pack / 26 packs
C1	1 yellow stripe	up to 75 cm	over 40.5 cm	up to 75 cm	23 580	1 pack / 26 packs
C2	2 yellow stripes	75-85 cm	over 40.5 cm	up to 75 cm	23 581	1 pack / 26 packs
C3	3 yellow stripes	over 85 cm	over 40.5 cm	up to 75 cm	23 582	1 pack / 26 packs

Knee length, one pair per pack

Size	Colour code foot	Knee length to hollow of the knee	Calf circumference	REF	Shipping Units
A	1 blue stripe	up to 40 cm	up to 34.5 cm	22 530	1 pack / 30 packs
B	1 green stripe	40-43 cm	34.5-40.5 cm	22 531	1 pack / 30 packs
C	1 yellow stripe	up to 47 cm	over 40.5 cm	22 532	1 pack / 30 packs