



# Durelast®

## Durelast® Combi very short stretch

Bandage with very short stretch, exceptionally high compressive power and durable elasticity; for strong compression of the extremities in phlebology as well

as for the treatment of acute and chronic lymphoedemas; for support and relief with distortions and contusions in traumatology and sports medicine.

### Durelast® Bandage

Firm-grip fabric, skin coloured, 66% cotton, 34% polyamide, with durable elasticity due to textured polyamide warp threads, very short stretch, extensibility approx. 45%, can be washed and boiled, sterilisable (A 121 °C).



<b>Durelast® Bandage</b>			
	Width	REF	Shipping Units
<b>5 m</b> (stretched length) with bandage clips individually in cardboard box	6 cm	22 100	1 pc. / 108 pcs.
	8 cm	22 101	1 pc. / 108 pcs.
	10 cm	22 102	1 pc. / 90 pcs.
	12 cm	22 103	1 pc. / 54 pcs.
<b>5 m</b> (stretched length) with bandage clips, unwrapped, 10 per box	6 cm	22 150	10 pcs. / 160 pcs.
	8 cm	22 151	10 pcs. / 80 pcs.
	10 cm	22 152	10 pcs. / 120 pcs.
	12 cm	22 153	10 pcs. / 80 pcs.
<b>Durelast® Combi</b>		Durelast bandages in varying widths in pairs, in cardboard box	
<b>5 m</b> (stretched length) with bandage clips	8 cm / 10 cm	22 110	1 box / 30 boxes