



Fortelast® • Fortelast® Plus

Fixona® Elastic Adhesive Bandages

Lengthwise elastic bandages, made of skin-friendly fabric, with excellent extensibility as a result of the porously applied

synthetic adhesive. Both bandages are permeable to air and moisture. They are indicated for compression and support bandaging in traumatology, func-

tional bandaging in sports medicine, and in phlebology for compression.



Fortelast®

White, 85% cotton, 15% viscose, lengthwise elastic, extensibility approx. 70%, hypoallergenic porously applied synthetic adhesive.

Fixona®

White, 100% cotton, lengthwise elastic, extensibility approx. 80%, synthetic adhesive with increased adhesive power.

| Fortelast® | Width | REF | Shipping Units |
|---|--------|--------|-----------------|
| 2.5 m (unstretched length) individually in cardboard box | 3 cm | 94 740 | 1 pc. / 75 pcs. |
| | 6 cm | 94 741 | 1 pc. / 90 pcs. |
| | 8 cm | 94 742 | 1 pc. / 90 pcs. |
| | 10 cm | 94 743 | 1 pc. / 72 pcs. |
| | 15 cm | 94 745 | 1 pc. / 24 pcs. |
| Fortelast® Plus | | | |
| rolled inwards | 2.5 cm | 16 752 | 1 pc. / 84 pcs. |
| | 5 cm | 16 753 | 1 pc. / 60 pcs. |
| | 7.5 cm | 16 754 | 1 pc. / 60 pcs. |
| | 10 cm | 16 755 | 1 pc. / 72 pcs. |
| | 15 cm | 16 756 | 1 pc. / 24 pcs. |
| rolled outwards | 2.5 cm | 16 757 | 1 pc. / 84 pcs. |
| | 5 cm | 16 758 | 1 pc. / 60 pcs. |
| | 7.5 cm | 16 759 | 1 pc. / 60 pcs. |
| | 10 cm | 16 760 | 1 pc. / 72 pcs. |
| Fixona® | | | |
| 4.5 m (stretched length) individually in cardboard box adhesive coating on the outside | 5 cm | 94 751 | 1 pc. / 50 pcs. |
| | 7.5 cm | 94 752 | 1 pc. / 50 pcs. |
| | 10 cm | 94 753 | 1 pc. / 50 pcs. |
| | 15 cm | 94 755 | 1 pc. / 50 pcs. |