

Idealbinde

short stretch

Bandage type “Ideal” conforming to DIN 61632, textile-elastic, the classic bandage for everyday use in hospital and doctor’s surgery; for powerful compression of the extremities in phlebology and lymphology; for pre-, intra-, and postoperative thrombo-prophylaxis; for sup-

port and relief in cases of distortions and contusions and for sport supports; for treating tendosynovitis, and for securing splints. “Ideal” bandages are air-permeable and skin-friendly; their firm-grip fabric structure ensures that the bandage turns do not slip.

Idealbinde

White fabric, conforming to DIN 61632, at least 85% cotton, viscose, textile-elastic through twisted cotton warp threads; short stretch, extensibility approx. 90%, can be washed and boiled.



Idealbinde			
	Width	REF	Shipping Units
5 m (stretched length) with bandage clips individually sealed	4 cm	14 100	1 pc. / 80 pcs.
	6 cm	14 101	1 pc. / 100 pcs.
	8 cm	14 102	1 pc. / 100 pcs.
	10 cm	14 103	1 pc. / 60 pcs.
	12 cm	14 104	1 pc. / 80 pcs.
	15 cm	14 105	1 pc. / 60 pcs.
	20 cm	14 106	1 pc. / 35 pcs.
	30 cm	14 108	1 pc. / 15 pcs.
5 m (stretched length) without bandage clips unwrapped 10 per box	4 cm	14 110	10 pcs. / 80 pcs.
	6 cm	14 111	10 pcs. / 100 pcs.
	8 cm	14 112	10 pcs. / 100 pcs.
	10 cm	14 113	10 pcs. / 60 pcs.
	12 cm	14 114	10 pcs. / 80 pcs.
	15 cm	14 115	10 pcs. / 60 pcs.
	20 cm	14 116	10 pcs. / 60 pcs.
	30 cm	14 118	10 pcs. / 40 pcs.