

In ready-for-use formats of different sizes, with lasting vertical elasticity; for compression in the treatment of oedemas in phlebology and lymphology; also used for functional bandages in sports medicine.

The fine-pored material with its velvety surface can also be easily cut to the required shape.

Free of natural latex.



## Komprex® Foam Rubber Pads

Synthetic latex (styrene-butadiene), permeable to air and water vapour.

## Komprex® Foam Rubber Bandages

Komprex bandages are indicated for severe oedemas, thrombophlebitis and after extensive sclerotherapy of varices. The bandages are applied in spiral turns, without stretching, and edge-to-edge – i.e. without overlapping; the increased compression effect is achieved by applying elastic bandages on top. If a

dressing is cut to a specific shape from a sheet of Komprex or from a bandage (e.g. U- and L-shaped pieces), the edges should be chamfered to prevent pressure spots.

Pads, bandages or material cut to size can also be used as partial padding with dressings and splints.

Synthetic latex (styrene-butadiene), permeable to air and water vapour.

## Komprex® II Foam Plate, wavy

Narrow strips of foam material, embedded in a skin-friendly nonwoven, offer an extremely breathable, non-slip padding material. The special, wavy surface structure of Komprex II creates a firm yet gentle micro-massaging effect.

<b>Komprex® Foam Rubber Pads</b>				
Size	Thickness	Size	REF	Shipping Units
individually in cardboard box				
00 oval	10 mm	9 x 6.5 cm	22 300	1 pc. / 60 pcs.
0 kidney-shaped	10 mm	9 x 5 cm	22 301	1 pc. / 75 pcs.
1 kidney-shaped	12 mm	12 x 6.5 cm	22 302	1 pc. / 30 pcs.
2 rectangular	12 mm	17 x 13 cm	22 303	1 pc. / 25 pcs.
4 trapezoid	5 mm	39/22 x 33 cm	22 305	1 pc. / 20 pcs.
<b>Komprex® Foam Rubber Bandages</b>				
	5 mm	8 cm x 1 m	22 310	1 pc. / 30 pcs.
rolled	10 mm	8 cm x 1 m	22 311	1 pc. / 20 pcs.
individually	10 mm	10 cm x 1 m	22 312	1 pc. / 20 pcs.
in cardboard box	5 mm	8 cm x 2 m	22 313	1 pc. / 30 pcs.
	10 mm	8 cm x 2 m	22 314	1 pc. / 15 pcs.
<b>Komprex® Foam Rubber Sheet</b>				
rolled, in plastic pouch	10 mm	100 x 50 cm	22 320	1 pc. / 10 pcs.
<b>Komprex® II Foam Plate, wavy</b>				
individually packed in a pouch		65 x 65 cm	16 985	1 pc. / 6 pcs.