

# Lenkelast®

## Elastra Ideal medium stretch

Firm-grip, all purpose bandages with durable elasticity and with a wide range of applications: everyday retention bandages with compressive wound care, or splint dressings, for support

and pressure relief; also for light compression.

The extensibility of these bandages makes them conform well to all body contours, and the dressing remains in posi-

tion. The fine fabric is air-permeable and highly skin-friendly.

### Lenkelast®

White, 64% cotton, 35% polyamide, 1% elastane, medium stretch, extensibility approx. 130%.

### Elastra Ideal

White, 58% cotton, 41% polyamide, 1% elastane, medium stretch, extensibility approx. 135%, washable, sterilisable (A 121°C).



| Lenkelast®   | Width | REF    | Shipping Units     |
|--|-------|--------|--------------------|
| 5 m<br>(stretched length)<br>individually sealed     | 6 cm  | 19 460 | 10 pcs. / 250 pcs. |
|  | 8 cm  | 19 461 | 10 pcs. / 200 pcs. |
|  | 10 cm | 19 462 | 10 pcs. / 240 pcs. |
|  | 12 cm | 19 463 | 10 pcs. / 200 pcs. |
|  | 15 cm | 19 464 | 10 pcs. / 160 pcs. |
|  | 20 cm | 19 465 | 10 pcs. / 100 pcs. |
| 5 m<br>(stretched length)<br>unwrapped<br>10 per box | 6 cm  | 19 490 | 10 pcs. / 250 pcs. |
|  | 8 cm  | 19 491 | 10 pcs. / 200 pcs. |
|  | 10 cm | 19 492 | 10 pcs. / 240 pcs. |
|  | 12 cm | 19 493 | 10 pcs. / 200 pcs. |
|  | 15 cm | 19 494 | 10 pcs. / 160 pcs. |
|  | 20 cm | 19 495 | 10 pcs. / 100 pcs. |
| Elastra Ideal  |       |        |                    |
| 5 m<br>(stretched length)<br>unwrapped<br>10 per box | 6 cm  | 10 122 | 10 pcs. / 160 pcs. |
|  | 8 cm  | 10 124 | 10 pcs. / 120 pcs. |
|  | 10 cm | 10 108 | 10 pcs. / 120 pcs. |
|  | 12 cm | 10 110 | 10 pcs. / 100 pcs. |
|  | 15 cm | 10 116 | 10 pcs. / 80 pcs.  |
|  | 20 cm | 10 118 | 10 pcs. / 60 pcs.  |