



# Lenkideal®

short stretch

Short stretch bandage with durable elasticity, very similar to the DIN "Ideal" bandage, but with longer-lasting elasticity and more durable, easy-care and economical – a bandage for all situations: for compressive

retention of wound dressings; for the application of ointment dressings; for compression of the extremities in phlebology, for pre-, intra- and postoperative thrombo-prophylaxis; for support and relief in cases of

distortions and contusions as well as for sport supports; for treating tendosynovitis, and for securing splints.

Lenkideal with its flat fabric does not become bulky, is air-permeable and skin-friendly.



## Lenkideal®

White, dense fabric, 64% cotton, 35% polyamide, 1% elastane, with durable elasticity, short stretch, extensibility approx. 90%.



Lenkideal®			
	Width	REF	Shipping Units
<b>5 m</b> (stretched length) with bandage clips individually sealed	6 cm	14 581	1 pc. / 100 pcs.
	8 cm	14 582	1 pc. / 140 pcs.
	10 cm	14 583	1 pc. / 120 pcs.
	12 cm	14 584	1 pc. / 120 pcs.
	15 cm	14 585	1 pc. / 120 pcs.
<b>5 m</b> (stretched length) without bandage clips unwrapped 10 per box	4 cm	19 580	10 pcs. / 140 pcs.
	6 cm	19 581	10 pcs. / 100 pcs.
	8 cm	19 582	10 pcs. / 140 pcs.
	10 cm	19 583	10 pcs. / 120 pcs.
	12 cm	19 584	10 pcs. / 120 pcs.
	15 cm	19 585	10 pcs. / 120 pcs.
	20 cm	19 586	10 pcs. / 60 pcs.