

Mollelast® haft

Mollelast® haft color

Cohesive conforming bandage with special stretch characteristics.

Mollelast haft is lengthwise elastic and widthwise extensible. This makes the bandage perfectly suitable for application on cone-shaped parts of the body, and offers a comfortable

and secure fit even in joint areas.

The new Mollelast haft was further optimised and provided with safe woven edges. The side threads are firmly anchored in the knit, so that they cannot become untied or fray.

Mollelast haft is also suitable as a simple support for sporting activities, as well as for cannula retention. If the dressing is fitted correctly, it will remain in position without cutting in, and is both air – permeable and well tolerated by the skin.



Mollelast® haft

Knitted bandage, white, 66% viscose, 34% polyamide, lengthwise elasticity approx. 65%, widthwise extensibility approx. 50%, made cohesive by a microfine latex coating.

Mollelast® haft color

Knitted bandage, red and blue, 56% viscose, 44% polyamide, lengthwise elasticity approx. 55%, widthwise extensibility approx. 50%, made cohesive by a microfine latex coating.

Mollelast® haft			
individually in cardboard box	Width	REF	Shipping Units
white 4 m (stretched length)	4 cm	30 063	1 pc. / 120 pcs.
	6 cm	30 064	1 pc. / 100 pcs.
	8 cm	30 065	1 pc. / 120 pcs.
	10 cm	30 066	1 pc. / 120 pcs.
	12 cm	30 067	1 pc. / 100 pcs.
white 20 m (stretched length)	4 cm	30 068	1 pc. / 18 pcs.
	6 cm	30 069	1 pc. / 18 pcs.
	8 cm	30 070	1 pc. / 18 pcs.
	10 cm	30 071	1 pc. / 18 pcs.
	12 cm	30 072	1 pc. / 18 pcs.
Mollelast® haft color			
individually in cardboard box			
red 20 m (stretched length)	6 cm	30 073	1 pc. / 18 pcs.
	8 cm	30 074	1 pc. / 18 pcs.
	10 cm	30 075	1 pc. / 18 pcs.
blue 20 m (stretched length)	6 cm	30 076	1 pc. / 18 pcs.
	8 cm	30 077	1 pc. / 18 pcs.
	10 cm	30 078	1 pc. / 18 pcs.