

# Perfekta® Cohesive

long stretch


Long stretch bandage with durable elasticity, cohesive on both sides; a support and relief bandage for the prophylactic relief of capsules and ligaments, prevention of haematomas, and for contusions, distortions, dislocations and tendon injuries; also used as a compressive bandage.

Its cohesivity makes it especially suitable for conically shaped or frequently moved parts of the body.



## Perfekta® Cohesive

Long stretch bandage with durable elasticity, skin coloured, 96% cotton, 4% elastane, lengthwise extensibility approx. 140%, cohesive on both sides due to microfine latex dots, washable at 40°C.



Perfekta® Cohesive			
individually in cardboard box	Width	REF	Shipping Units
<b>5 m</b> (stretched length)	6 cm	20 106	1 pc. / 50 pcs.
	8 cm	20 107	1 pc. / 50 pcs.
	10 cm	20 104	1 pc. / 50 pcs.
	12 cm	20 105	1 pc. / 50 pcs.