

Porelast® • Porelast Superforte

Porelast® PrO₂ Adhesive Bandages

Lengthwise elastic adhesive bandages, with a choice of two types of adhesive, for permanent compression bandages in the field of phlebology and for functional bandages in sports medicine; for support and relief with distortions, contusions and

dislocations; also for redressment bandages, e.g. for malposition of the foot.

The fabric with porously applied adhesive is permeable to air and water vapour, thus offering good ventilation even with overlapping bandage turns.



Porelast® Porelast® Superforte

Elastic textile fabric, with selvedges, 100% cotton, skin coloured, lengthwise extensibility approx. 80%, strong synthetic rubber adhesive, generally well tolerated by the skin, X-ray-permeable.

The difference between the two bandages is the higher adhesive power of the Superforte version.

Porelast® PrO₂

Elastic textile fabric with cut edges, 100% cotton, skin coloured, lengthwise extensibility approx. 70%, hypoallergenic adhesive, with high air-permeability, X-ray-permeable.

Porelast®	Width	REF	Shipping Units
2.5 m (unstretched length) individually in cardboard box	6 cm	30 210	1 pc. / 50 pcs.
	8 cm	30 211	1 pc. / 50 pcs.
	10 cm	30 212	1 pc. / 50 pcs.
2.5 m (unstretched length) unwrapped	6 cm	30 213	10 pcs. / 50 pcs.
	8 cm	30 214	10 pcs. / 50 pcs.
	10 cm	30 215	10 pcs. / 50 pcs.
Porelast® Superforte			
2.5 m (unstretched length) individually in cardboard box	6 cm	95 410	1 pc. / 50 pcs.
	8 cm	95 411	1 pc. / 50 pcs.
	10 cm	95 412	1 pc. / 50 pcs.
Porelast® PrO₂			
2.5 m (unstretched length) individually in cardboard box	6 cm	30 220	1 pc. / 50 pcs.
	8 cm	30 221	1 pc. / 50 pcs.
	10 cm	30 222	1 pc. / 50 pcs.
2.5 m (unstretched length) unwrapped	6 cm	30 223	10 pcs. / 50 pcs.
	8 cm	30 224	10 pcs. / 50 pcs.
	10 cm	30 225	10 pcs. / 50 pcs.