

Porodress®

Adhesive Bandage

Widthwise elastic adhesive bandage, lengthwise inelastic, with strong adhesive power.

Porodress adhesive bandages are used chiefly as a stirrup for lesions on capsules or ligaments and for muscle and tendon injuries or disorders; also for redressement and traction,

especially in combination with Porelast applied in circular turns.

The fabric with porously applied adhesive is permeable to air and water vapour, thus offering good ventilation even with overlapping bandage turns.



Porodress®

Elastic textile fabric, with cut edges, skin coloured, 58% cotton, 42% viscose, widthwise extensibility approx. 60%, strong synthetic rubber adhesive, generally well tolerated by the skin, X-ray-permeable.

Porodress®			
individually in cardboard box	Width	REF	Shipping Units
2.5 m (lengthwise inelastic)	6 cm	30 240	1 pc. / 50 pcs.
	8 cm	30 241	1 pc. / 50 pcs.
	10 cm	30 242	1 pc. / 50 pcs.