

Porotape® Haftan®

Porotape®

Non-elastic adhesive tape with high tensile strength, which can be torn both lengthwise and widthwise, with an effective zinc oxide rubber adhesive.

Porotape can be applied to the skin directly, but is generally used in combination with adhesive bandages, for partial immobilization with a support and

relieving function on joints, muscles and ligaments; also as a prophylactic dressing aimed at avoiding overloading.

The adhesive tape is dispensed easily from the reel; the serrated edge makes tearing easy and prevents frayed edges.

Fine fabric, 100% viscose, zinc oxide rubber adhesive, raw white.



Haftan®

Extremely thin bandage (0.8 mm) made from polyurethane foam, open-pored, supple. Used as skin protection under adhesive bandages and tape dressings, as well as for securing pads which provide extra compression. The material's surface makes bandage turns adhere well to each other.

Haftan is also used for the smooth application of padding material underneath P.o.P. and synthetic casts.

Polyurethane foam, white, resistant to ointments, sweat and cleaning agents.

Porotape®		REF	Shipping Units
length: 10 m	2.5 cm	30 081	6 pcs. / 120 pcs.
	3.8 cm	30 084	1 pc. / 96 pcs.
	3.8 cm	30 082	12 pcs. / 96 pcs.
	5.0 cm	30 083	6 pcs. / 72 pcs.
Haftan®			
thickness: 0.8 mm	7 cm x 27.5 m	22 315	6 pcs. / 60 pcs.
	10 cm x 27.5 m	22 316	6 pcs. / 60 pcs.