

Rosidal® K • Elko®-Rosidal® K

Raucodur Kräftig (strong) short stretch

Short stretch bandages with high compressive power and durable elasticity; for strong compression of the extremities in phlebology as well as for the treatment of acute and chronic

lymphoedemas; for support and relief in traumatology and sports medicine.

The above bandages are permeable to air and well tolerated by the skin. The firm-grip fabric

texture ensures that the bandage turns hold securely.

The cohesive short stretch bandage Rosidal haft (see page 81) offers even greater protection against slipping.



Rosidal® K Raucodur Kräftig (strong)

Firm-grip fabric, skin coloured, 100% cotton, textile-elastic due to twisted warp threads, short stretch, extensibility approx. 90%, washable up to 95°C.

Elko®-Rosidal® K

For compressive dressings using the double bandaging technique, two Rosidal bandages of 10 cm width are available under the name of Elko-Rosidal K.

Rosidal® K			
	Width	REF	Shipping Units
5 m (stretched length) with bandage clips (except the bandage with 4 cm) individually in cardboard box	4 cm	22 199	1 pc. / 125 pcs.
	6 cm	22 200	1 pc. / 108 pcs.
	8 cm	22 201	1 pc. / 108 pcs.
	10 cm	22 202	1 pc. / 90 pcs.
	12 cm	22 203	1 pc. / 54 pcs.
10 m (stretched length) without bandage clips individually in cardboard box	10 cm	22 206	1 pc. / 60 pcs.
	12 cm	22 205	1 pc. / 40 pcs.
5 m (stretched length) with bandage clips unwrapped, 10 per box	8 cm	22 250	10 pcs. / 140 pcs.
	10 cm	22 251	10 pcs. / 140 pcs.
	12 cm	22 252	10 pcs. / 120 pcs.
Elko®-Rosidal® K 2 Rosidal K bandages packed in pairs in cardboard box			
5 m (stretched length) with bandage clips	10 cm	22 225	1 box / 30 boxes
Raucodur Kräftig (strong)			
5 m (stretched length) with bandage clips individually in cardboard box	6 cm	20 168	1 pc. / 50 pcs.
	8 cm	20 169	1 pc. / 50 pcs.
	10 cm	20 164	1 pc. / 50 pcs.
	12 cm	20 166	1 pc. / 50 pcs.
unwrapped 10 per box	8 cm	20 170	10 pcs. / 100 pcs.
	10 cm	20 165	10 pcs. / 100 pcs.
	12 cm	20 167	10 pcs. / 100 pcs.
10 m	10 cm	20 268	30 pcs.