

Application	Lymphology Set 1 for leg, small REF 59 795	Lymphology Set 2 for leg, large REF 59 796	Lymphology Set 3 for arm, small REF 59 797	Lymphology Set 4 for arm, large REF 59 798	Lymphology Set 5 for leg, small REF 16 838	Lymphology Set 6 for leg, large REF 16 839
tg® Tubular Bandage For skin protection (hygiene) under decongestive compression dressings.	Size 7 (1x)	Size 9 (1x)	Size 5 (1x)	Size 7 (1x)	Size 9 (1x)	Size 9 (1x)
Mollelast® Bandage Comfortably soft, highly elastic conforming bandage for the compression of fingers and toes in lymphological compression dressings.	4 cm x 4 m (20x)	4 cm x 4 m (20x)	4 cm x 4 m (20x)	4 cm x 4 m (20x)	4 cm x 4 m (20x)	6 cm x 4 m (20x)
Komprex® Foam Rubber Bandage To cut out U- and L-shaped pieces (edges should be chamfered) for padding and increasing pressure in the ankle joint. Additionally, as basic material for padding packs (strips or balls) and for the local softening of tissue.	10 cm x 1 m Thickness 1 cm (1x)	10 cm x 1 m Thickness 1 cm (1x)	10 cm x 1 m Thickness 1 cm (1x)	10 cm x 1 m Thickness 1 cm (1x)	8 cm x 1 m Thickness 0.5 cm (1x)	10 cm x 1 m Thickness 1 cm (1x)
Cellona® Synthetic Padding Padding under compression bandages for the even distribution of pressure and for the prevention of constrictions and pressure spots.	6 cm x 3 m (8x) 10 cm x 3 m (8x)	10 cm x 3 m (8x) 15 cm x 3 m (8x)	6 cm x 3 m (4x) 10 cm x 3 m (4x)	6 cm x 3 m (4x) 10 cm x 3 m (4x)		
Rosidal® soft Skin-friendly, air-permeable foam bandage for padding under compression dressings. Prevents the dressing from slipping.					10 cm x 2.5 m Thickness 0.3 cm (2x) 15 cm x 2.5 m Thickness 0.4 cm (2x)	15 cm x 2.5 m Thickness 0.3 cm (2x) 15 cm x 2.5 m Thickness 0.4 cm (3x)
Rosidal® K Bandage (strong) Short stretch bandage, 100% cotton, for strong compression of extremities in lymphological treatment.	6 cm x 5 m (2x) 8 cm x 5 m (2x) 10 cm x 5 m (2x) 12 cm x 5 m (4x)	6 cm x 5 m (2x) 8 cm x 5 m (2x) 10 cm x 5 m (4x) 12 cm x 5 m (6x)	6 cm x 5 m (2x) 10 cm x 5 m (2x) 12 cm x 5 m (2x)	6 cm x 5 m (2x) 10 cm x 5 m (2x) 12 cm x 5 m (3x)	6 cm x 5 m (2x) 8 cm x 5 m (2x) 10 cm x 5 m (2x) 12 cm x 5 m (4x)	6 cm x 5 m (2x) 8 cm x 5 m (2x) 10 cm x 5 m (4x) 12 cm x 5 m (6x)
Idealbinde For strong compression and dressing stabilisation on the thigh and abdomen.					20 cm x 5 m (1x)	20 cm x 5 m (2x)