

Rosidal® haft short stretch

Raucodur Cohesive

Cohesive short stretch bandages with strong compression; for compression of the extremities with disorders of the veins and lymphatic system, as permanent bandages, for the pre- and postoperative thromboprophylaxis, and for compression of the thigh to complement the lower leg bandage. The compression remains therapeutically effective over a prolonged period.

These bandages are also suitable for use as functional bandages, for immobilising, supporting and relieving capsules, ligaments, muscles and joints.

The stretch characteristics are similar to those of Rosidal K (see page 89), and – despite their adhesive properties – the bandages are permeable to air and water vapour.



Skin coloured fabric, 100% cotton, elastic textile with twisted warp threads, low stretch, cohesive through microfine latex dots; can be washed at 40°C.

Rosidal haft: Extensibility approx. 60%.

Raucodur Cohesive: Extensibility approx. 75%.



| Rosidal® haft | | | |
|----------------------------------|-------|--------|-----------------|
| individually in cardboard box | Width | REF | Shipping Units |
| 5 m (stretched length) | 6 cm | 22 270 | 1 pc. / 60 pcs. |
| | 8 cm | 22 271 | 1 pc. / 40 pcs. |
| | 10 cm | 22 272 | 1 pc. / 40 pcs. |
| | 12 cm | 22 273 | 1 pc. / 30 pcs. |
| Raucodur Cohesive | | | |
| individually in cardboard box | | | |
| 5 m (stretched length) | 8 cm | 20 173 | 1 pc. / 50 pcs. |
| | 10 cm | 20 171 | 1 pc. / 50 pcs. |
| | 12 cm | 20 172 | 1 pc. / 48 pcs. |