

Rosidal® soft

Rosidal soft is ideally suited for padding underneath lymphological and phlebological compression bandages. The padding bandage leads to an even pressure distribution of short stretch bandages (e.g. Rosidal K) over the extremity, this makes the compression treatment safe by preventing constrictions. The bandage will not slip as the open pores of Rosidal soft are interlocking.

The skin-friendly, latex-free bandage made from polyurethane foam is highly permeable to air and moisture.

Rosidal soft is particularly tear-resistant and does not lose its elasticity even after repeated washing.

(We recommend a washing temperature of 40 - 60°C. If the bandage is extremely dirty, it may be boiled.)



100% polyurethane foam, highly extensible, washable up to 95°C, sterilisable (A 121°C, A 134°C).



Rosidal® soft			
	Size	REF	Shipping Units
in bulk (unwrapped, with bandage clips)	10 cm x 0.3 cm x 2.5 m	23 100	1 box / 36 pcs.
	10 cm x 0.4 cm x 2.5 m	23 101	24 pcs.
	12 cm x 0.4 cm x 2.5 m	23 102	18 pcs.
	15 cm x 0.4 cm x 2.5 m	23 103	14 pcs.
individually packed (in cardboard box, with bandage clips)	10 cm x 0.3 cm x 2.5 m	23 110	1 pc. / 30 pcs.
	10 cm x 0.4 cm x 2.5 m	23 111	1 pc. / 24 pcs.
	12 cm x 0.4 cm x 2.5 m	23 112	1 pc. / 16 pcs.
	15 cm x 0.4 cm x 2.5 m	23 113	1 pc. / 14 pcs.
two pieces in a cardboard box (with bandage clips)	10 cm x 0.2 cm x 2 m	13 016	1 box / 23 boxes