

Rosidal® sys

Short Stretch Compression System

For an effective compression therapy of venous ulcers.

The use of this set in long-term treatment is particularly effective and economical. All components are carefully adjusted to

each other, are of high and lasting quality and are combined in a handy set.

The reusability of the padding and compression bandages results in a considerable cost-ef-

fectiveness compared with disposable products. The padding bandage Rosidal soft and the short stretch bandage Rosidal K are washable up to 50 times at up to 95°C.

A system solution also saves handling time and offers highest therapy safety. The administration of individual components is eliminated.

For wound treatment, the dressing can be adequately added, so that the wound is optimally cared for.

It is possible to schedule bandage changes individually for each patient. Thus, a compression bandage can be worn for up to 7 days.



Rosidal® sys			
Components	Size		REF
Set	Shipping Units 1 piece/4 pieces		16 984
All components may be ordered individually (except those marked with *).			
tg® tubular bandage For skin protection under the decongesting compression dressing	Size 7	1x	24 005
Rosidal® soft bandage Skin-friendly, latex-free foam bandage made from 100% polyurethane foam; washable up to 95°C	10 cm x 2 m x 0.2 cm	4x	13 016
Rosidal® K bandage Tightly compressing short stretch bandage made from 100% cotton; washable up to 95°C	8 cm x 5 m	2x	22 201
	10 cm x 5 m	2x	22 202
Porofix® Skin coloured adhesive tape on the roll; coated with synthetic rubber adhesive	2.5 cm x 5 m	1x	34 311
Mollelast® haft bandage Soft cohesive conforming bandage, lengthwise and widthwise elastic	8 cm x 20 m	1x	30 070
Laundry bag*		1x	