

# Varicex® Zinc Paste Bandages (Unna's Boot) Starched Gauze Bandages

Bandages with different fabric texture, impregnated with zinc paste; ready-for-use, time-saving application of semi-rigid permanent bandages for the treatment of thrombophlebitis, oedemas, and chronic venous insufficiency. Also used for leg ulcers in the healing phase and for the

aftercare of fractures. Whether they are applied in circular turns or with the semi-figure-of-eight technique – Varicex zinc paste bandages do not bulk up.

## Varicex® S

Bi-elastic knitted bandage for application in circular turns. It is par-

ticularly easy to apply, because of the material's stretch effect.

## Varicex® E

Elastic creped bandage for application in circular turns.

## Varicex® F

Produced from gauze bandages, with a low moisture content, for rapid drying.

## Varicex® T

Produced from gauze bandages, with a high moisture content, for slow drying.

## Starched Gauze Bandage

20-thread, starched, for light rigid bandaging or for the reinforcement of bandages with an immobilising function.



<b>Varicex® S</b>		
individually in aluminium pouch	REF	Shipping Units
10 cm x 5 m (stretched length)	20 630	1 pc. / 80 pcs.
10 cm x 7 m (stretched length)	20 631	1 pc. / 60 pcs.
10 cm x 10 m (stretched length)	20 632	1 pc. / 30 pcs.
individually in aluminium pouch, 20 pcs. in a cardboard box		
10 cm x 5 m (stretched length)	20 635	20 pcs.
10 cm x 7 m (stretched length)	20 636	20 pcs.
10 cm x 10 m (stretched length)	20 637	20 pcs.
<b>Varicex® E</b> individually in aluminium pouch		
8 cm x 5 m (stretched length)	20 620	1 pc. / 60 pcs.
<b>Varicex® F</b> individually in aluminium pouch		
10 cm x 5 m	20 600	1 pc. / 80 pcs.
10 cm x 7 m	20 601	1 pc. / 60 pcs.
10 cm x 10 m	20 602	1 pc. / 50 pcs.
<b>Varicex® T</b> individually in aluminium pouch		
10 cm x 10 m	20 612	1 pc. / 50 pcs.
<b>Starched Gauze Bandage</b> unwrapped, 10 per box		
6 cm x 4 m	14 701	10 pcs. / 100 pcs.
8 cm x 4 m	14 702	10 pcs. / 100 pcs.
10 cm x 4 m	14 703	10 pcs. / 100 pcs.
12 cm x 4 m	14 704	10 pcs. / 50 pcs.
15 cm x 4 m	14 705	10 pcs. / 50 pcs.