



tg[®] grip Elasticated Tubular Support Bandage

Tubular bandage with durable elasticity, supple, exerting even pressure; highly practical and handy for dressing retention with light compression, for support and relief on ligaments,

tendons and joints and for wide-area (thorax) dressings.

The bandage is air-permeable, slightly padding and very skin-friendly through its high cotton

content. tg grip can be fitted, removed and re-applied easily and without creasing. The compression effect can be increased with a second layer or by selecting a smaller size.



tg[®] grip Elasticated Tubular Support Bandage

Seamlessly knitted tube, raw white, 85% cotton, 10% elastodiene (rubber), covered with 5% polyamide yarn, widthwise elastic, can be washed (60°C) and sterilised (A 134°, EO, irradiation).



tg[®] grip Elasticated Tubular Support Bandage

10 m roll in cardboard dispenser	Size	REF	Shipping Units
for children's hands, children's forearms	A (4.6 cm)	24 120	1 roll / 24 rolls
for children's arms, children's lower legs	B (6.0 cm)	24 121	1 roll / 16 rolls
for slim arms, children's legs	C (6.5 cm)	24 122	1 roll / 16 rolls
for medium-sized arms, lower legs	D (7.5 cm)	24 123	1 roll / 16 rolls
for large arms, legs, slim thighs	E (8.75 cm)	24 124	1 roll / 12 rolls
for legs, medium-sized thighs	F (10 cm)	24 125	1 roll / 12 rolls
for legs, for large thighs	G (12 cm)	24 126	1 roll / 8 rolls
for small trunks	J (17.5 cm)	24 127	1 roll / 6 rolls
for medium-sized trunks	K (21.5 cm)	24 128	1 roll / 4 rolls
for large trunks	L (32.5 cm)	24 129	1 roll / 4 rolls
for extra-large trunks	M (37.5 cm)	24 130	1 roll / 2 rolls