

Cellacare® Cervical

Cellacare® Cervical Plus

Anatomically shaped cervical collars, alternatively with or without reinforcement, for the posttraumatic and postoperative treatment of the cervical spine, as well as for the relief of the neck muscles.

The hook & loop closure makes the cervical collars easy to handle and infinitely adjustable. The removable protective cotton cover can be washed at 60°C,

making the collars suitable for repeated use over an extended period of time.

Cellacare® Cervical Collar

Protective cover, stockinette cover, core and hook & loop closure identical with Cellacare Cervical Plus; shape: anatomically shaped curve in the chin area, slightly higher in the neck.

Cellacare® Cervical Plus Orthosis

Protective cover: cotton, polyamide, elastodien; stockinette cover: double-layer cotton; core: polyurethane foam; hook & loop closure: polyamide; reinforcement: polyethylene; shape: anatomically shaped curve in the chin area, slightly higher in the neck.

Care Instructions:
protective cover: machine washable up to 60°C, cervical collars: hand wash recommended, do not iron, do not tumble-dry, not suitable for dry-cleaning. Colour: ecru.



Cellacare® Cervical individually in a cardboard box		
Height	REF	Shipping Units
6 cm	22 490	1 pc. / 20 pcs.
8 cm	22 491	1 pc. / 16 pcs.
10 cm	22 492	1 pc. / 14 pcs.
Cellacare® Cervical Plus individually in a cardboard box		
6 cm	22 495	1 pc. / 20 pcs.
8 cm	22 496	1 pc. / 16 pcs.
10 cm	22 497	1 pc. / 14 pcs.