

Cellacare® Cervidur Schanz design

Wrap-around cervical collar available in three different widths, tried and trusted for support and stabilisation in the region of the cervical vertebral column. The cervical collar relieves pressure on the neck muscles, promotes recovery. Due to its immobilizing and

warming effect it aids medicinal treatment and physiotherapy in cases of cervical vertebral column syndrome, plexus neuralgia and torticollis.

Handling is simple. Cervidur cervical collars are skin-friendly and easy-care.

Cellacare® Cervidur

Firm, stable nonwoven core (polyester, polypropylene) with seamless knitted stockinette (cotton), with longitudinal stitching and ties (cotton); off-white; can be washed at up to 40°C, does not require ironing.



Cervidur® Cervical Collar

rolled, individually in a cardboard box

Size	REF	Shipping Units
6 cm x 1.80 m	22 405	1 pc. / 40 pcs.
8 cm x 1.80 m	22 406	1 pc. / 32 pcs.
10 cm x 1.80 m	22 407	1 pc. / 24 pcs.