

Cellacare® Clavicula

Ready-to-wear support on the rucksack brace principle; for the conservative treatment of fractures of the clavicle. Also post-operatively for immobilisation of the shoulder region with extension of the clavicle.

Three buckles and flat hook & loop closures make putting on and adjustment fast and simple.

A soft pad cushions the back strap and buckle junction area.

The skin-friendly materials of this support provide a high level of wearer comfort.

Cellacare® Clavicula

Belts: foam padding: polyurethane, outer cover and quilted seams: cotton, hook & loop tapes: polyamide,

velour strips: cotton, buckles: polyamide;

Back padding: needle-punched nonwoven wool, viscose, knitted cover sheathing and seam material: cotton.

Washable up to 30°C, do not iron, do not tumble-dry, not suitable for dry-cleaning, colour: white.



Cellacare® Clavicula

individually packed in a cardboard box

Size	Shoulder Width	REF	Shipping Units
0	18 to 25 cm	22 420	1 pc. / 54 pcs.
1	25 to 33 cm	22 421	1 pc. / 54 pcs.
2	33 to 40 cm	22 422	1 pc. / 54 pcs.
3	40 to 50 cm	22 423	1 pc. / 54 pcs.
4	50 to 60 cm	22 424	1 pc. / 54 pcs.