

# Cellacare® Dorsal

Back support with elastic construction elements and an unobtrusive optical appearance for the moderate stabilization and relief of the lumbar region. Suitable for the treatment and prevention of acute and chronic lumbalgia and lumbal ischialgia, in cases of false posture and irritation, as well as in the early stage of osteoporosis.

The elastic support provides gentle massage to the back and abdomen. Flexible stays in the back ensure the required stabi-

lization. This proprioceptively acting support stimulates the metabolism, reduces muscle tension and assists the neuromuscular activity.

The high quality materials are air permeable and prevent the accumulation of moisture or heat. The trapezoidal spacer fabric in the back offers temperature balancing and a comfortable padding and support function.

The compressive action and thereby the degree of stabilization can easily be adjusted to individual needs.



## Cellacare® Dorsal

Materials: cotton, polyester, elastodiene, polyamide, polypropylene.

Colour: white/grey

Washable up to 30°C, do not tumble-dry, iron or dry-clean, no chlorine bleach.



## Cellacare® Dorsal

individually packed in a cardboard box

Size	Waist Circumference	REF	Shipping Units
1	60 - 75 cm	17 992	1 pc. / 10 pcs.
2	75 - 95 cm	17 993	1 pc. / 10 pcs.
3	95 - 115 cm	17 994	1 pc. / 12 pcs.
4	115 - 135 cm	17 996	1 pc. / 12 pcs.
5	> 135 cm	17 997	1 pc. / 12 pcs.