

Cellacare® Dotop F

Cellacare® Dotop M



Back support with stabilising stays in the back and a special pulley mechanism (EasyFit) for stabilisation and stress relief in the lumbar region; available in 2 versions, for the male and the female anatomy. The product comes with additional exercise bands for optional therapy to strengthen the muscles of the upper body.

The new EasyFit mechanism is based on the “pulley” principle - in closing the straps, an optimum support effect is achieved with a minimum of effort on the patient’s part.

With its special reinforcement in the lumbar region, Cellacare®

Dotop prevents extreme misalignment (lordosis/kyphosis) of the lumbar spine. Functional mobilisation of the patient can be facilitated by removing the outer stabilising stays in the back.

The new Dotop therapy concept is consistently in compliance with the current recommendation to treat back problems not by rest, but by specific movements and exercises which improve the functionality of the back muscles. Moreover, this approach can efficiently reduce the level of pain.

Cellacare Dotop is custom-tailored not only to the specific

male/female anatomical differences, it also takes into account the different aesthetic preferences of male and female patients. The two versions are available in 5 sizes each. The support is very easy to apply.

Indications:

- therapy and prevention of acute and chronic lumbar syndrome (lumbalgia, lumboschialgia, protrusion, prolapse, facet syndrome, spondylolysis)
- chronic poor posture and dysfunction

Cellacare® Dotop

Support: cotton, polyester, elastodiene, polyamide, polypropylene, elastane, viscose.

Exercise bands: polypropylene, polyester, polyamide.

Washable at temperatures of up to 30°C; do not tumble-dry; do not iron; not suitable for dry-cleaning; no chlorine bleaching.

Colour: Cellacare Dotop F: black with pink; Cellacare Dotop M: grey/black with orange.



Cellacare® Dotop F



Cellacare® Dotop M

Cellacare® Dotop F

Back support for women

Size	Waist Circumference	REF	Shipping Units
1	65 - 75 cm	25 541	1 pc. / 16 pcs.
2	75 - 90 cm	25 542	1 pc. / 16 pcs.
3	90 - 105 cm	25 543	1 pc. / 16 pcs.
4	105 - 120 cm	25 544	1 pc. / 16 pcs.
5	120 - 135 cm	25 545	1 pc. / 16 pcs.

Cellacare® Dotop M

Back support for men

1	65 - 75 cm	25 546	1 pc. / 16 pcs.
2	75 - 90 cm	25 547	1 pc. / 16 pcs.
3	90 - 105 cm	25 548	1 pc. / 16 pcs.
4	105 - 120 cm	25 549	1 pc. / 16 pcs.
5	120 - 135 cm	25 550	1 pc. / 16 pcs.