

Cellacare® Epi

Active elbow support with inserts and a stabilising strap, for stress relief in tendon insertions on the elbow; for the therapy of lateral and medial epicondylitis, postoperative and posttraumatic irritation and excessive stress on the elbow joint, as well as effusion into the joint and swelling caused by arthrosis or arthritis.

Firm on the outside and soft on the inside - due to the specific material composition and three-dimensional knitting.

The three-fold therapy system with its highly compressive active knit, integrated muscle-relieving inserts, and circular stabilising strap has a massaging and mobilising effect which promotes blood circulation and reduces swelling.

The stabilising strap features an in-woven stretch indicator for additional increase of pressure and individual adjustment as required.

The excellent fit, soft edges, and the special soft comfort zone in the crook of the arm provide superior wearer comfort



Cellacare® Epi

Suitable for right and left arm, three-dimensional knitting, cotton/nylon fabric, inserts: Komplex, insert pockets: cotton, stabilising strap.

Washable at temperatures of up to 30 °C, do not tumble-dry, do not iron, not suitable for dry cleaning. Colour: blue/black.

Cellacare® Epi			
Size	Arm Circumference (9 cm from elbow to distal)	REF	Shipping Units
0	18-20 cm	13 361	1 pc. / 20 pcs.
1	20-22 cm	13 362	1 pc. / 20 pcs.
2	22-24 cm	13 363	1 pc. / 20 pcs.
3	24-26 cm	13 364	1 pc. / 20 pcs.
4	26-28 cm	13 365	1 pc. / 20 pcs.
5	28-30 cm	13 366	1 pc. / 20 pcs.
6	30-32 cm	13 367	1 pc. / 20 pcs.