

Cellacare® Genu

Active knee support for the functional treatment of posttraumatic and postoperative conditions of the knee; patella femoral pain syndrome, myopathy/tendopathy, ligament instability, effusion into the joint, and swelling caused by active arthrosis or arthritis. Also in cases of contusion or ligament damage, and for injury prevention.

Firm on the outside and soft on the inside - due to the specific material composition and three-dimensional knitting.

Anatomically shaped at the knee joint, with diagonal stretch towards the top and bottom edges.

The integrated ring insert (made of silicone) provides an intermittent massaging effect and improved physiological distribution of pressure.

Cellacare Genu provides support and compression of the knee joint. It promotes blood circulation, stimulates the metabolism, and increases mobility, thus reducing effusion and swelling, and relieving pain.



Highly breathable

NEW:

- improved widthwise and lengthwise elasticity on the thigh
 - less slipping of the support
 - improved fit and wearer comfort
- the seam in the back is offset to the side
- popliteal insert
 - improved fit

Cellacare® Genu

Three-dimensional knitting, colour: blue/black, cotton/nylon fabric, insert: silicone, insert pocket: cotton

Washable at temperatures of up to 30 °C, do not tumble-dry, do not iron, not suitable for dry cleaning

Cellacare® Genu

individually packed in a cardboard box

Size	Circumference 12 cm below center of knee	Circumference 14 cm above center of knee	REF	Shipping Units
0	25 - 28 cm	35 - 38 cm	16 849	1 pc. / 36 pcs.
1	28 - 31 cm	38 - 41 cm	16 850	1 pc. / 36 pcs.
2	31 - 34 cm	41 - 44 cm	16 851	1 pc. / 36 pcs.
3	34 - 37 cm	44 - 47 cm	16 852	1 pc. / 36 pcs.
4	37 - 40 cm	47 - 50 cm	16 853	1 pc. / 36 pcs.
5	40 - 43 cm	50 - 53 cm	16 854	1 pc. / 36 pcs.
6	43 - 46 cm	53 - 56 cm	16 855	1 pc. / 36 pcs.