

# Cellacare® Genucast

Sleeve-shaped knee splint for immobilisation of the joint in cases of inflammation and knee injuries of various types, also postoperatively in a functional, slightly bent position, individually adjustable.

Aluminum back and side stays (bent 15°), which can be removed from their pockets, secure and stabilize. The hook &

loop closed retaining straps around the leg are quick and easy to use.

The skin-friendly cotton lining is padded with open-cell foam padding. Cellacare Genucast, the “gentle” knee splint, combines the securing effect of a rigid brace with wearer comfort and the hygienic advantages of a removable textile orthosis.

## Cellacare® Genucast

Foam padding: polyurethane/polyester; white inside lining: cotton, blue outside cover: polyester, edging: cotton, hook & loop tapes and eyelets: polyamide, longitudinal tapes: polypropylene, sewing thread: cotton polyester, aluminum stays.

Washable up to 30°C, remove stays before washing, do not iron, do not tumble-dry, not suitable for dry-cleaning, colour: navy/grey.



### Cellacare® Genucast

individually in PE bag

Size	for Body Size	Height	Thigh Circumf.	Calf Circumf.	REF	Shipping Units
S	< 1.45 m	40 cm	30 to 38 cm	26 to 31 cm	22 445	1 pc. / 17 pcs.
M	1.45 to 1.55 m	48 cm	35 to 43 cm	28 to 32 cm	22 446	1 pc. / 14 pcs.
L	1.56 to 1.75 m	54 cm	44 to 54 cm	32 to 38 cm	22 447	1 pc. / 14 pcs.
XL	1.76 to 1.85 m	60 cm	53 to 62 cm	39 to 44 cm	22 448	1 pc. / 29 pcs.
XXL	> 1.85 m	66 cm	56 to 68 cm	42 to 48 cm	22 449	1 pc. / 26 pcs.